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INSTRUCTIONS FOR USE

Kit for specimen collection, transport, and RNA isolation

Bi-CoV[®] set

Intended Use

Bi-CoV[®] Set is a sampling kit used for a non-invasive collection of biological material from the front nose or oral cavity by soaking saliva into the inserted swab and subsequent specimen transport to the laboratory for PCR diagnostics. Alternatively, nasopharyngeal or oropharyngeal swab specimens can be used as biological material. The set is suitable for self-sampling. The collection and transport solution ensures a gradual decomposition of the viral particles, release of the nucleic acids (RNA) into the solution and their protection against degradation. **The sample is ready for PCR analysis without the need to perform RNA isolation in the diagnostic laboratory.**


Packaging:

10 sets (Cat. No. BI001-10) 25 sets (Cat. No. BI001-25L)
 25 sets (Cat. No. BI001-25)
 50 sets (Cat. No. BI001-50)

For one sampling:

- 2 ml plastic vial (BI001-25L: 10–13 ml plastic test tube) with 1 ml of collection/transport solution
- swab for sampling saliva or from the nasal walls of the nostrils (or nasopharyngeal / oropharyngeal) sampling

Instructions for sampling from the front nose or saliva collection from the oral cavity*:

Collection from the front nose:

1. Before sampling, the nose should be blown.
2. Open the package (tear the foil and paper on the opposite side from the sampling swab) to remove the plastic rod with the swab.
3. Open the plastic tube with the solution by turning the cap.
4. Insert the swab 2 cm deep while rubbing gently the nasal walls of the nostrils 3x in each direction.

Collection from the oral cavity (saliva):

1. Tested person must not brush teeth, use mouthwash, drink, eat, use chewing gum or smoke at least 30 minutes before the sample collection.
2. The tested person sniffs and clears the throat to get mucus and saliva from the nasopharynx and throat gets to the oral cavity.
3. Open the package (tear the foil and paper on the opposite side from the swab) to remove the plastic rod with the swab.
4. Open the plastic tube with the solution by turning the cap.
5. Insert the swab into the oral cavity to the collected saliva place (on or under the tongue, near the molars, etc.), wipe around the palatoglossal arches).

- Insert the swab with the collected biological material into the tube with the solution, break the plastic rod and close the cap tightly.
- Label the tube and send it to the diagnostic laboratory
- The sampling must be performed carefully to obtain a representative and analysable material. While collecting the sample, use personal protective equipment. Direct contact with the biological material can be dangerous.

*From **nasopharynx or oropharynx**, the collection is performed by a circular movement of the swab on the posterior wall of the nasopharynx or oropharynx.

Storage and other information:

- Storage possible at room temperature if used up within 2 months, in other cases store at 2-8 °C. Avoid exposure to direct sunlight.
- For immediate use after opening.
- The product can be used as long as the packaging is not damaged.
- For single use only.
- Dispose of the used material as infectious waste.

Instructions for sample preparation for PCR test:

- The sample is ready for PCR testing upon receipt in the laboratory.
- The required amount of sample can be taken without having to remove the swab.
- To increase the sensitivity of the PCR test, a short centrifugation is recommended (10,000-12,000 g, 3-5 minutes)
- The sample is compatible with majority of the PCR kits, recommended:

Manufacturer	Cat. No.	Name
GeneProof a.s.	COV2/GP/100	GeneProof SARS-CoV-2 PCR Kit
	COV2S/GP/100	GeneProof SARS-CoV-2 Screening PCR Kit
ELISABETH PHARMACON, spol. s r.o.	90077-RT-A	EliGene [®] COVID19 BASIC A RT
BAG Diagnostics GmbH	728250	ViroQ SARS-CoV-2
GENERI BIOTECH s.r.o.	3231-200	gb SARS-CoV-2 Multiplex
iNtRON Biotechnology (dodavatel Novatin)	1000056 ER	Lilif [™] COVID-19 Real-time RT-PCR Kit
Shanghai ZJ Bio-Tech	RR-0479-02	Liferiver Novel Coronavirus (2019-nCoV) Real Time Singleplex

Manufacturer: Bioinova, a.s.
 Videnska 1083
 142 00 Prague 4, Czech Republic
www.bioinova.cz

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Bi-CoV® set

Intended Use

Bi-CoV® Set is a sampling kit used for a non-invasive collection of biological material from the front nose or oral cavity by soaking saliva into the inserted swab and subsequent specimen transport to the laboratory for PCR diagnostics. Alternatively, nasopharyngeal or oropharyngeal swab specimens can be used as biological material. The set is suitable for self-sampling. The collection and transport solution ensures a gradual decomposition of the viral particles, release of the nucleic acids (RNA) into the solution and their protection against degradation. **The sample is ready for PCR analysis without the need to perform RNA isolation in the diagnostic laboratory.**

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
***From nasopharynx or oropharynx**, the collection is performed by a circular movement of the swab on the posterior wall of the nasopharynx or oropharynx.

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