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## INSTRUCTIONS FOR USE

### Kit for specimen collection, transport, and RNA isolation

### Bi-CoV®

#### Intended Use

Bi-CoV® is a sampling kit used for a non-invasive collection of biological material from the front nose or oral cavity by soaking saliva into the inserted swab and subsequent specimen transport to the laboratory for PCR diagnostics. Alternatively, nasopharyngeal or oropharyngeal swab specimens can be used as biological material. The set is suitable for self-sampling. The collection and transport solution ensures a gradual decomposition of the viral particles, release of the nucleic acids (RNA) into the solution and their protection against degradation. **The sample is ready for PCR analysis without the need to perform RNA isolation in the diagnostic laboratory.**


#### Packaging:

10 sets (Cat. No. BI001-10GP)  
50 sets (Cat. No. BI001-50GP)

#### For one sampling:

- 2 ml plastic vial (BI001-25L: 10–13 ml plastic test tube) with 1 ml of collection/transport solution
- swab for sampling saliva or from the nasal walls of the nostrils (or nasopharyngeal / oropharyngeal) sampling

#### Instructions for sampling from the front nose or saliva collection from the oral cavity\*:

#### Collection from the front nose:

1. Before sampling, the nose should be blown.
2. Open the package (tear the foil and paper on the opposite side from the sampling swab) to remove the plastic rod with the swab.
3. Open the plastic tube with the solution by turning the cap.
4. Insert the swab 2 cm deep while rubbing gently the nasal walls of the nostrils 3x in each direction.

#### Collection from the oral cavity (saliva):

1. Tested person must not brush teeth, use mouthwash, drink, eat, use chewing gum or smoke at least 30 minutes before the sample collection.
2. The tested person sniffs and clears the throat to get mucus and saliva form the nasopharynx and throat gets to the oral cavity.
3. Open the package (tear the foil and paper on the opposite side from the swab) to remove the plastic rod with the swab.
4. Open the plastic tube with the solution by turning the cap.
5. Insert the swab into the oral cavity to the collected saliva place (on or under the tongue, near the molars, etc.), wipe around the palatoglossal arches).

- Insert the swab with the collected biological material into the tube with the solution, break the plastic rod. # If necessary, the collected material on the swab can be rinsed into the collection solution by a circular movement in the solution for at least 10 seconds. The swab can then be removed and discarded. Close the cap tightly.
- Label the tube and send it to the diagnostic laboratory
- The sampling must be performed carefully to obtain a representative and analysable material.
- ATTENTION: Direct contact with the biological material can be dangerous. While collecting the sample, use personal protective equipment.

\*From **nasopharynx** or **oropharynx**, the collection is performed by a circular movement of the swab on the posterior wall of the nasopharynx or oropharynx.

#### Storage and other information:

- Storage possible at room temperature if used up within 2 months, in other cases store at 2-8 °C. Avoid exposure to direct sunlight.
- For immediate use after opening.
- The product can be used as long as the packaging is not damaged.
- For single use only.
- Dispose of the used material as infectious waste.

#### Instructions for sample preparation for PCR test:

- The sample is ready for PCR testing upon receipt in the laboratory.
- The required amount of sample can be taken without having to remove the swab.
- To increase the sensitivity of the PCR test, a short centrifugation is recommended (10,000-12,000 g, 3-5 minutes).
- Serious incidents possibly related to the product (e.g. serious deterioration of health) shall be reported to the manufacturer and national competent authority.

Manufacturer	Cat. No.	Name
GeneProof a.s.	COV2/GP/100	GeneProof SARS-CoV-2 PCR Kit
	COV2S/GP/100	GeneProof SARS-CoV-2 Screening PCR Kit
	COV2A/GP/100	GeneProof SARS-CoV-2 Advanced PCR Kit

**Manufacturer:** Bioinova, a.s.,  
Videnska 1083  
142 00 Prague 4, Czech Republic  
[www.bioinova.cz](http://www.bioinova.cz)