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Literatura:

1. **Zav'yalov VP et al.** 2019. Interferon-Inducible Myxovirus Resistance Proteins: Potential Biomarkers for Differentiating Viral from Bacterial Infections. Clin Chem. 65(6):739-750.
2. **Michell PS et al.** 2014. An evolutionary perspective on the broad antiviral specificity of MxA. Curr Opin Microbiol. 16(7):493–499.
3. **Toivonen L et al.** 2015. Blood MxA protein as a marker for respiratory virus infections in young children. J Clin Vir. 62: 8-13.
4. **Schneider-Schaulies S et al.** 1994. Cell type-specific MxA-mediated inhibition of measles virus transcription in human brain cells. J Virol. 68(11):6910-6917.
5. **Frese M et al.** 1996. Inhibition of bunyaviruses, phleboviruses, and hantaviruses by human MxA protein. J Virol. 70(2):915-923.
6. **Nakabayashi M et al.** 2006. MxA-based recognition of viral illness in febrile children by a whole blood assay. Pediatr Res. 60(6):770-4.
7. **Ronni T et al.** 1993. Control of IFN-inducible MxA gene expression in human cells. J Immunol. 150(5):1715-1726.
8. **Abrahams ME et al.** 1995. IFN-α induces MxA gene expression in cultured human corneal fibroblasts. Exp Eye Res. 60(2):137-142.
9. **Engelmann I et al.** 2015. Diagnosis of viral infections using myxovirus resistance protein A (MxA). Pediatrics. 135:e985–93.
10. **Sambursky R et Shapiro N.** 2015. Evaluation of a combined MxA and CRP point-of-care immunoassay to identify viral and/or bacterial immune response in patients with acute febrile respiratory infection. Eur Clin Respir J. 2:28245.

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Literatur:

- Zav'yalov VP et al.** 2019. Interferon-Inducible Myxovirus Resistance Proteins: Potential Biomarkers for Differentiating Viral from Bacterial Infections. Clin Chem. 65(6):739-750.
- Michell PS et al.** 2014. An evolutionary perspective on the broad antiviral specificity of MxA. Curr Opin Microbiol. 16(7):493–499.
- Toivonen L et al.** 2015. Blood MxA protein as a marker for respiratory virus infections in young children. J Clin Vir. 62: 8-13.
- Schneider-Schaulies S et al.** 1994. Cell type-specific MxA-mediated inhibition of measles virus transcription in human brain cells. J Virol. 68(11):6910-6917.
- Frese M et al.** 1996. Inhibition of bunyaviruses, phleboviruses, and hantaviruses by human MxA protein. J Virol. 70(2):915-923.
- Nakabayashi M et al.** 2006. MxA-based recognition of viral illness in febrile children by a whole blood assay. Pediatr Res. 60(6):770-4.
- Ronni T et al.** 1993. Control of IFN-inducible MxA gene expression in human cells. J Immunol. 150(5):1715-1726.
- Abrahams ME et al.** 1995. IFN-α induces MxA gene expression in cultured human corneal fibroblasts. Exp Eye Res. 60(2):137-142.
- Engelmann I et al.** 2015. Diagnosis of viral infections using myxovirus resistance protein A (MxA). Pediatrics. 135:e985–93.
- Sambursky R et Shapiro N.** 2015. Evaluation of a combined MxA and CRP point-of-care immunoassay to identify viral and/or bacterial immune response in patients with acute febrile respiratory infection. Eur Clin Respir J. 2:28245.

Hersteller: Bioinova, a.s.
Videnska 1083
142 00 Prag 4 – Krc
Tschechische Republik
www.bioinova.cz

INSTRUCTIONS FOR USE
Rapid test for the detection of viral infections
Bi-VirTest®

Intended use:
Bi-VirTest® is a rapid quantitative immunochromatographic test that helps to distinguish the cause of an incipient respiratory infection. The test is intended for near-patient testing.

Bi-VirTest® measures the level of MxA (myxovirus resistance protein A) in a capillary blood sample. MxA is produced in cells in response to elevated levels of antiviral interferons type I and III (1). Among other mediators of the immune response, it is significant for its activity against an unusually broad spectrum of viruses (2). Elevated MxA levels have been observed during infections with rhinovirus, influenza virus, human parainfluenza virus, coronavirus, respiratory syncytial virus, or human metapneumovirus (3). Its antiviral activity also includes some viruses causing infections beyond the respiratory tract (4; 5).

Clinical diagnosis is based on the fact that normal MxA level in blood is very low and increases specifically in response to acute viral infection (6). MxA increase occurs within hours and the highest concentration is reached within the first day of infectious disease (7; 8). In the presence of interferon, elevated levels are maintained. MxA is also very stable (half-life ~ 2.3 days). Bacterial infections do not increase MxA level (6). The clinical utility of this biomarker in the determination of infections of viral origin has been confirmed by numerous clinical studies (3; 9; 6; 10).

Indications:
The test is intended for patients with clinical symptoms of acute respiratory infection. It should be performed when symptoms are present, at least 6–12 hours after the onset of fever and preferably within the first three days of its onset.


Packaging:
10 tests (catalog number BI005-10)

- Contents:**
- | | | | |
|--------------------------------|-----|------------------------------|----|
| - Test cassette | 10x | - Wash Buffer in Unisampler™ | 1x |
| - UniSampler™ | 10x | - RFID card | 1x |
| - Lysis Buffer in 3–5mL bottle | 1x | - Instructions for use | 1x |


 **WARNING:** Contact with buffers (Lysis Buffer and Wash Buffer) may cause severe eye irritation. Contact for emergency situations (UK customers): National Poisons Information Service, Royal Infirmary of Edinburgh, Edinburgh EH16 4SA, tel.: +44 121 507 4123, 844 892 0111.

- Materials and equipment required but not provided:**
- Lancets 23G, gloves, sterile swabs, and disinfectant for blood collection
 - Reader: Bi-Reader® (opTricon GmbH, DE) or Igloo Reader Pro (Dx365 GmbH, DE)


Test procedure:


UniSampler™ contains 2 parts: container and cap with a double capillary for blood collection.


Add 6 drops of **Lysis Buffer** into the container of the empty UniSampler™.


Obtain a drop of capillary blood using the standard procedure. Touch the tip of the capillary to the blood droplet.


Wait until both parts of the capillary fill with blood.


Insert the full capillary into the container and press firmly to close the container completely.


Shake the UniSampler™ by vigorous vertical motion to move the blood to the container and mix it with Buffer.


Remove the cap.

Invert UniSampler™ and gently squeeze 2 drops onto the sample well on the test cassette.

Wait for 4 minutes and then apply 2 drops of **Wash Buffer** to the sample well.

Wait for 11 minutes (incubation time) and read the result using one of the readers.

Bi-Reader®


Igloo Reader Pro

 **WARNING:** Wash Buffer must be applied exactly 4 to 5 minutes after the sample, otherwise the result may be invalid.






The result-reading procedure for **Bi-Reader®** is provided in the section below.
The result-reading procedure for **Igloo Reader Pro** is provided on page 9.

Bi-Reader® procedure:

Place the cassette adapter onto the used test cassette. The pin on the bottom of the adapter must fit into the sample well.

Make sure that the detection window on the bottom of the reader is clean and place the reader on the cassette adapter. The slanted corner of the reader must face its counterpart on the adapter.

1. Turn on the reader by pressing the button on the front. The reader runs a self-test, during which "WAIT" is displayed. If the reader can display the last result, it will now appear on the display. Confirm the result by <u>briefly</u> (< 1 s) pressing the button. After the beep, "ON" is displayed. Press the button again <u>briefly</u> *.
 → 
2. "RFID" appears on the display. Place a lot-specific RFID card provided with the test onto the top of the reader. This will upload the calibration data to the reader.
 → 
3. After the beep signal, "TEST" is displayed. Press the button to display "RUN". After the next beep signal, the MxA concentration in ng/mL is displayed.
 → 

*Note: If you want to use the automatic readout function after the incubation time of the test (11 minutes), press the button longer (>1 s) in step 1. After step 3, a countdown indication appears on the display (i.e., a repeated countdown from 9 to 0) and the result is displayed automatically after the incubation time. If you press the button during the countdown, the measurement is terminated.

Interpretation of results:

Test result (ng/mL)	MxA level	Acute viral infection	Symptoms of the disease	Recommendations about further diagnosis and use of antibiotics*
0–45	Normal (negative)	NO	Mild	Primarily NO, then according to the course of the disease + CRP test may be done.
			Severe	YES (CRP test immediately + use of antibiotics)
>45	Elevated (positive)	YES	Mild	Primarily NO, then according to the course of the disease + CRP test may be done.
			Severe	

*** MxA is an important marker of active viral infection. However, when making the final diagnosis, it is essential to consider the patient's clinical condition, results of any additional tests, and the possibility of a mixed infection, which may cause elevated levels of both MxA and CRP.**

Limitations:

The test may not be accurate in the following circumstances:

- Fever lasting more than 7 days.
- Immunodeficiency, or use of immunosuppressive drugs.
- Use of oral medications for infection.
- Within 30 days after vaccination with live virus vaccines.
- Use of interferons (e.g., to treat multiple sclerosis).
- Within 30 days after a major injury or surgery.
- Presence of atypical bacterial respiratory infections (e.g. *Mycoplasma pneumoniae* or *Chlamydophila pneumoniae*).
- Presence of active infections outside the respiratory tract (e.g. herpesviruses).

Analytical performance:

Sensitivity

The analytical sensitivity (lower limit of detection) is 20 ng/mL.

Measurement range

The measurement range is 45 ng/mL to 2 000 ng/mL. The calibration curve contained 13 points and the samples were measured in 10-fold determination.

Interfering substances

Interference was measured on MxA positive and negative samples in 5-fold determination. Low interference indicates a bias of more than 20%, medium interference more than 30% and high interference more than 50%.

Name of substance	Concentration	Interference	Name of substance	Concentration	Interference
Ampicillin-Na	1 mg/mL	None	Cyclosporine	5 µg/mL	Low
Bilirubin, conjugated	0.6 mg/mL	Low	Doxycycline Hyclate	5 µg/mL	Low
Bilirubin, unconjugated	0.66 mg/mL	None	Levodopa	30 µg/mL	None
Biotin	3.5 µg/mL	High	Phenylbutazone	0.1 mg/mL	None
	1.2 µg/mL	Medium	Rifampicin	60 µg/mL	None
	0.6 µg/mL	Low	Theophylline	0.1 mg/mL	None
	0.1 µg/mL	Low			

Clinical performance:

Sensitivity and specificity

In a clinical study, 43 patients with symptoms of acute upper respiratory tract infection were evaluated. MxA levels obtained by Bi-VirTest® were compared with the evidence of viral infection in a respiratory sample by PCR (reference test). The resulting sensitivity and specificity are 91.7% and 94.7%, respectively.

		Reference test		
		Positive	Negative	Total
Bi-VirTest®	Positive	22	1	23
	Negative	2	18	20
	Total	24	19	43

Storage and other information:

- Appropriate protective equipment (gloves) must be worn when working with the test and during blood collection.
- Store at 2-8 °C, out of direct sunlight.
- Keep out of the reach of children.
- The cassette must be used immediately after opening.
- Open buffers (Lysis Buffer and Wash Buffer) must be stored under the storage conditions specified for the entire kit and can be used for the entire shelf life.
- The product can be used only if the packaging is not damaged. Do not use damaged parts of the test.
- The test cassette and empty UniSampler™ are intended for single use.
- Used test cassettes must be disposed of as infectious waste. Residual buffers must be disposed of as hazardous waste.
- Read the result only with the Bi-Reader® and its accessories recommended by the test manufacturer.
- Any serious adverse events related to the product (e.g., serious deterioration of health) shall be reported to the manufacturer and local competent authority if required by local regulations.

Literature:

1. **Zav'yalov VP et al.** 2019. Interferon-Inducible Myxovirus Resistance Proteins: Potential Biomarkers for Differentiating Viral from Bacterial Infections. Clin Chem. 65(6): 739-750.
2. **Michell PS et al.** 2014. An evolutionary perspective on the broad antiviral specificity of MxA. Curr Opin Microbiol. 16(7):493–499.
3. **Toivonen L et al.** 2015. Blood MxA protein as a marker for respiratory virus infections in young children. J Clin Vir. 62: 8-13.
4. **Schneider-Schaulies S et al.** 1994. Cell type-specific MxA-mediated inhibition of measles virus transcription in human brain cells. J Virol. 68(11):6910-6917.
5. **Frese M et al.** 1996. Inhibition of bunyaviruses, phleboviruses, and hantaviruses by human MxA protein. J Virol. 70(2):915-923.
6. **Nakabayashi M et al.** 2006. MxA-based recognition of viral illness in febrile children by a whole blood assay. Pediatr Res. 60(6):770-4.
7. **Ronni T et al.** 1993. Control of IFN-inducible MxA gene expression in human cells. J Immunol. 150(5):1715-1726.
8. **Abrahams ME et al.** 1995. IFN-α induces MxA gene expression in cultured human corneal fibroblasts. Exp Eye Res. 60(2):137-142.
9. **Engelmann I et al.** 2015. Diagnosis of viral infections using myxovirus resistance protein A (MxA). Pediatrics. 135:e985–93.
10. **Sambursky R et Shapiro N.** 2015. Evaluation of a combined MxA and CRP point-of-care immunoassay to identify viral and/or bacterial immune response in patients with acute febrile respiratory infection. Eur Clin Respir J. 2:28245.

Producer: Bioinova, a.s.
Videnska 1083
142 00 Prague 4 – Krc
Czech Republic
www.bioinova.cz

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
Litteratur:

1. **Zav'yalov VP et al.** 2019. Interferon-Inducible Myxovirus Resistance Proteins: Potential Biomarkers for Differentiating Viral from Bacterial Infections. Clin Chem. 65(6):739-750.
2. **Michell PS et al.** 2014. An evolutionary perspective on the broad antiviral specificity of MxA. Curr Opin Microbiol. 16(7):493-499.
3. **Toivonen L et al.** 2015. Blood MxA protein as a marker for respiratory virus infections in young children. J Clin Vir. 62: 8-13.
4. **Schneider-Schaulies S et al.** 1994. Cell type-specific MxA-mediated inhibition of measles virus transcription in human brain cells. J Virol. 68(11):6910-6917.
5. **Frese M et al.** 1996. Inhibition of bunyaviruses, phleboviruses, and hantaviruses by human MxA protein. J Virol. 70(2):915-923.
6. **Nakabayashi M et al.** 2006. MxA-based recognition of viral illness in febrile children by a whole blood assay. Pediatr Res. 60(6):770-4.
7. **Ronni T et al.** 1993. Control of IFN-inducible MxA gene expression in human cells. J Immunol. 150(5):1715-1726.
8. **Abrahams ME et al.** 1995. IFN-α induces MxA gene expression in cultured human corneal fibroblasts. Exp Eye Res. 60(2):137-142.
9. **Engelmann I et al.** 2015. Diagnosis of viral infections using myxovirus resistance protein A (MxA). Pediatrics. 135:e985-93.
10. **Sambursky R et Shapiro N.** 2015. Evaluation of a combined MxA and CRP point-of-care immunoassay to identify viral and/or bacterial immune response in patients with acute febrile respiratory infection. Eur Clin Respir J. 2:28245.

Produsent: Bioinova, a.s.
Videnska 1083
142 00 Praha 4 - Krc
Tsjeckia
www.bioinova.cz

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Bibliografia:

1. **Zav'yalov VP et al.** 2019. Interferon-Inducible Myxovirus Resistance Proteins: Potential Biomarkers for Differentiating Viral from Bacterial Infections. Clin Chem. 65(6):739-750.

2. **Michell PS et al.** 2014. An evolutionary perspective on the broad antiviral specificity of MxA. Curr Opin Microbiol. 16(7):493–499.

3. **Toivonen L et al.** 2015. Blood MxA protein as a marker for respiratory virus infections in young children. J Clin Vir. 62: 8-13.

4. **Schneider-Schaulies S et al.** 1994. Cell type-specific MxA-mediated inhibition of measles virus transcription in human brain cells. J Virol. 68(11):6910-6917.

5. **Frese M et al.** 1996. Inhibition of bunyaviruses, phleboviruses, and hantaviruses by human MxA protein. J Virol. 70(2):915-923.

6. **Nakabayashi M et al.** 2006. MxA-based recognition of viral illness in febrile children by a whole blood assay. Pediatr Res. 60(6):770-4.

7. **Ronni T et al.** 1993. Control of IFN-inducible MxA gene expression in human cells. J Immunol. 150(5):1715-1726.

8. **Abrahams ME et al.** 1995. IFN-α induces MxA gene expression in cultured human corneal fibroblasts. Exp Eye Res. 60(2):137-142.

9. **Engelmann I et al.** 2015. Diagnosis of viral infections using myxovirus resistance protein A (MxA). Pediatrics. 135:e985–93.

10. **Sambursky R et Shapiro N.** 2015. Evaluation of a combined MxA and CRP point-of-care immunoassay to identify viral and/or bacterial immune response in patients with acute febrile respiratory infection. Eur Clin Respir J. 2:28245.

Producent: Bioinova, a.s.
Videnska 1083
142 00 Praha 4 – Krc
Republika Czeska
www.bioinova.cz

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EN

Igloo Reader Pro Procedure

1. Turn on the reader by pressing the power button located on the circular rubber bottom of the device.
2. Select “Start Measurement”. You can now enter the patient ID. Continue by pressing the ✓ symbol. Then you can enter the patient’s age and gender or deactivate the timer function*. Continue by pressing “Next”.
3. Insert the test cassette into the adapter supplied with the reader. Ensure the correct orientation of the cassette according to the markings on the adapter.
4. Insert the adapter with the test cassette into the reader.
5. After the incubation time (11 minutes) has elapsed, the measurement will start automatically. Do not remove the adapter from the reader during incubation or measurement.
6. Once the measurement is completed, the test result will be displayed on the screen.

*If the timer function is deactivated, the incubation time (11 minutes) must be measured separately. The measurement starts when the adapter with the cassette is inserted into the reader.

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